


















If the instructions are not followed properly, operations may lead to hazards for the product or the user/patient.

1.2. On the product/packaging

 Serial number	 Catalogue number
 Manufacturer	 Date of manufacture
 Class II (AC Adapter)	 Consult the Instructions for use
 This shows the Type B applied part.	 Do not dispose it with normal household waste
 Recycling	IPX0 Ordinary equipment
 EU-representative	 Atmospheric pressure for storage
 Temperature limits	 Storage humidity range
 CE marking	 Keep dry